

TC02_L17000_v1.31 - Agilent VEE Pro

File Edit View Debug Flow Device System I/O Data Display Tools Database Excel Window Help

Program Explorer

- TC02_L17000_v1.31.vee
 - Main
 - Initialize
 - RunTC02Sample
 - Inject
 - PumpAcid
 - Rinse_GetSample
 - Local UserFunctions
 - AcidFlush
 - Close_Dump
 - ErrorRecord
 - Full_Flow
 - Get_LicorData
 - Normal_Flow
 - Open_Dump
 - Plot
 - SampleSetUp
 - SaveData
 - Stop_Flow
 - Globals

Properties

Inject (UserObject) - M.215.196

- Appearance
 - Title: Inject
- Behavior
 - Breakpoint: False
 - TrigMode: Degrees
- Design
 - Declared Name
 - Declared Scope: Local to Context
 - ObjectIndex: M.215.196
- Popup Panel
 - PopupMoveable: True
 - PopupPanelTitle: UserObject
 - PopupPanelTitle: ActiveCaption
 - PopupPanelTitle Object Title Text
 - PopupPanelTitle: ActiveCaptionTex
 - ShowPanelOnE: False
 - ShowPopupPan: True
 - ShowPopupPan: True
- Title
 - Title string of the object

Inject

Note Pad

Standard is in position 8; rotate valve clockwise and inject.
Rotate counterclockwise and inject for any other sample.

SampleType

Name Data

SampleType

If/Then/Else

InjectStandard

/1o1V1000a0m0R

SyringeCom

InjectSample

/1o-1V1000a0m0R

SyringeCom

Serial10 (@ 10)

WRITE TEXT a EOL
< Double-Click to Add Transaction >

Error

Call ErrorRecord

ExecMode: VEE6&7 PROF MOD WEB