

TCO2\_L7000\_v1.31 - Agilent VEE Pro

File Edit View Debug Flow Device System I/O Data Display Tools Database Excel Window Help

Program Explorer

- TCO2\_L7000\_v1.31.vee
  - Main
    - Initialize
    - RunTCO2Sample
      - Inject
      - PumpAcid
      - Rinse\_GetSample
    - Local UserFunctions
      - AcidFlush
      - Close\_Dump
      - ErrorRecord
      - Full\_Flow
      - Get\_LicorData
      - Normal\_Flow
      - Open\_Dump
      - Plot
      - SampleSetUp
      - SaveData
      - Stop\_Flow
    - Globals

Properties

Serial12 (@ 12) (Direct I/O) - F21.0

**Appearance**

- BackColor ☐ Object
- Font ☐ Object Text
- ForeColor ☒ Object Text
- ShowTerminals True
- ShowTitleBar True
- Title Serial12 (@ 12)
- Title BackColor ☐ Direct I/O
- Title Font ☐ Object Title Text
- Title ForeColor ☒ Object Title Text

**Behavior**

- Breakpoint False

**Design**

- Declared Name
- Declared Scope Local to Context
- ObjectIndex F21.0

**Icon Appearance**

- IconPicture "to icon", Centered
- ShowIconTitle True

**Title**

Title string of the object

Serial12 (@ 12)

```
WRITE BYTE 254
WRITE BYTE 2
WRITE BYTE 254
WRITE BYTE 3
< Double-Click to Add Transaction >
```

Note Pad

```
Turn off relay #3
Turn off relay #4
```

Watch Terminals

Ready

ExecMode: VEE6&7 PROF MOD WEB